

Bone Screws – Demo

Cortical, cancellous, cannulated and headless compression screws.

Headless Compression Screw

A dual-pitch headless screw that buries beneath the articular surface and generates continuous compression, ideal for scaphoid and small-fragment fractures.

[Read online](#)

Cortical Bone Screw

Fully threaded self-tapping cortical screws for diaphyseal fixation and plate attachment in dense cortical bone.

[Read online](#)

Cancellous Bone Screw

Coarse-thread cancellous screws with deep threads for maximum purchase in metaphyseal and epiphyseal cancellous bone.

[Read online](#)

Cannulated Cancellous Screw

Guidewire-cannulated cancellous screws for percutaneous, accurately guided fixation of femoral neck and periarticular fractures.

[Read online](#)

Contact

- **Address:** 14 Harley Medical Court, London W1G 9PF, United Kingdom
- **Telephone:** [+91 9319821922](tel:+919319821922)
- **Email:** exports@orthomedimplants.com