

# Intramedullary Nails – Demo

Femoral, tibial, humeral and pediatric IM nails.

## Antegrade Femoral Interlocking Nail

A cannulated, closed-section titanium nail for antegrade fixation of femoral shaft fractures, engineered with an anatomic mediolateral bow and anti-rotation locking for rotational stability.

[Read online](#)

## Retrograde Femoral Nail

Designed for retrograde insertion through the intercondylar notch, ideal for distal femoral, periprosthetic and polytrauma fractures where antegrade access is unsuitable.

[Read online](#)

## Tibia Interlocking Nail

A proximal-bend tibial nail with multiplanar proximal locking for the full spectrum of tibial shaft and metaphyseal fractures, including suprapatellar insertion.

[Read online](#)

## Humerus Interlocking Nail

A low-profile humeral nail for proximal, shaft and combined humeral fractures, with a rounded proximal end that minimises rotator-cuff irritation.

[Read online](#)

## Proximal Femoral Nail (PFN)

A two-screw cephalomedullary nail for pertrochanteric and subtrochanteric femoral fractures, providing controlled compression across the fracture with anti-rotation support.

[Read online](#)

## Proximal Femoral Nail Anti-Rotation (PFNA)

A single helical-blade cephalomedullary nail that compacts cancellous bone for superior purchase in osteoporotic pertrochanteric fractures.

[Read online](#)

## Gamma-Type Trochanteric Nail

A short trochanteric nail with a large-diameter lag screw for stable and unstable intertrochanteric fractures, allowing early weight bearing.

[Read online](#)

## Expert Suprapatellar Tibial Nail

Optimised for semi-extended suprapatellar insertion, extending intramedullary fixation to very proximal and very distal tibial fractures.

[Read online](#)

## Reconstruction Femoral Nail

A reconstruction nail with two cephalic screws directed into the femoral neck for ipsilateral neck-and-shaft fractures and pathological lesions.

[Read online](#)

## Contact

- **Address:** 14 Harley Medical Court, London W1G 9PF, United Kingdom
- **Telephone:** [+91 9319821922](tel:+919319821922)
- **Email:** [exports@orthomedimplants.com](mailto:exports@orthomedimplants.com)