

Plates & Fixators – Demo

Locking compression plates, DHS/DCS, anatomical plates.

Dynamic Hip Screw (DHS) Plate

The gold-standard sliding hip screw for stable intertrochanteric and femoral neck fractures, allowing controlled fracture impaction under load.

[Read online](#)

Dynamic Condylar Screw (DCS) Plate

A fixed 95° condylar screw-and-plate construct for distal femoral and subtrochanteric fractures requiring rigid angular stability.

[Read online](#)

Proximal Humerus Locking Plate (PHILOS)

An anatomically contoured proximal humerus plate with a divergent locking-screw array and suture holes for two-, three- and four-part fractures.

[Read online](#)

Distal Femur Locking Plate (LISS)

A lateral distal femur plate designed for minimally invasive submuscular insertion (MIPO), acting as an internal fixator for comminuted supracondylar fractures.

[Read online](#)

Proximal Tibia Locking Plate

An anatomic lateral tibial plateau plate supporting the articular surface with a raft of subchondral locking screws for split-depression fractures.

[Read online](#)

Distal Radius Volar Locking Plate

A low-profile volar plate for distal radius fractures with a variable-angle distal row that lets the surgeon aim subchondral pegs beneath the articular surface.

[Read online](#)

Reconstruction Plate

A notched plate that can be bent and twisted in three planes for pelvic, acetabular, clavicle and periarticular fractures with complex geometry.

[Read online](#)

One-Third Tubular Plate

A thin, low-profile plate for lateral malleolus, olecranon and other subcutaneous locations where soft-tissue coverage is limited.

[Read online](#)

Clavicle Locking Plate

A superior, anatomically pre-contoured clavicle plate for displaced midshaft and lateral clavicle fractures, reducing the need for intra-operative bending.

[Read online](#)

Calcaneal Locking Plate

A perimeter calcaneal plate that supports the posterior facet and tuberosity in displaced intra-articular calcaneal fractures via a lateral extensile or sinus-tarsi approach.

[Read online](#)

Contact

- **Address:** 14 Harley Medical Court, London W1G 9PF, United Kingdom
- **Telephone:** [+91 9319821922](tel:+919319821922)
- **Email:** exports@orthomedimplants.com