

Implant Products – Demo

Headless Compression Screw

A dual-pitch headless screw that buries beneath the articular surface and generates continuous compression, ideal for scaphoid and small-fragment fractures.

[Read online](#)

Cannulated Cancellous Screw

Guidewire-cannulated cancellous screws for percutaneous, accurately guided fixation of femoral neck and periarticular fractures.

[Read online](#)

Cancellous Bone Screw

Coarse-thread cancellous screws with deep threads for maximum purchase in metaphyseal and epiphyseal cancellous bone.

[Read online](#)

Cortical Bone Screw

Fully threaded self-tapping cortical screws for diaphyseal fixation and plate attachment in dense cortical bone.

[Read online](#)

Monolateral External Fixator

A rigid modular monolateral frame for temporary damage-control fixation and definitive treatment of open tibial and femoral fractures.

[Read online](#)

Calcaneal Locking Plate

A perimeter calcaneal plate that supports the posterior facet and tuberosity in displaced intra-articular calcaneal fractures via a lateral extensile or sinus-tarsi approach.

[Read online](#)

Clavicle Locking Plate

A superior, anatomically pre-contoured clavicle plate for displaced midshaft and lateral clavicle fractures, reducing the need for intra-operative bending.

[Read online](#)

One-Third Tubular Plate

A thin, low-profile plate for lateral malleolus, olecranon and other subcutaneous locations where soft-tissue coverage is limited.

[Read online](#)

Reconstruction Plate

A notched plate that can be bent and twisted in three planes for pelvic, acetabular, clavicle and periarticular fractures with complex geometry.

[Read online](#)

Distal Radius Volar Locking Plate

A low-profile volar plate for distal radius fractures with a variable-angle distal row that lets the surgeon aim subchondral pegs beneath the articular surface.

[Read online](#)

Contact

- **Address:** 14 Harley Medical Court, London W1G 9PF, United Kingdom
- **Telephone:** [+91 9319821922](tel:+919319821922)
- **Email:** exports@orthomedimplants.com