


PLATES & FIXATORS

Calcaneal Locking Plate

Product code: CALC-PL



Specifications

Specification	Value
Material	Titanium Ti-6Al-4V
Diameter range	3.5mm locking screws
Length range	50mm – 90mm
Thread type	Perimeter locking
Availability	Available on order

Available sizes and configurations

- Small
- Medium
- Large

Product description

A perimeter calcaneal plate that supports the posterior facet and tuberosity in displaced intra-articular calcaneal fractures via a lateral extensile or sinus-tarsi approach.

Key Features

- Anatomic lateral wall contour
- Peripheral + posterior facet screws
- Fits extensile and sinus-tarsi approaches
- Very low profile to reduce wound issues

Technical Specifications

Product Code	CALC-PL
Material	Titanium Ti-6Al-4V
Diameter	3.5mm locking screws
Length Range	50mm – 90mm
Available Sizes	Small, Medium, Large
Thread / Locking	Perimeter locking
Certifications	ISO 13485, CE Mark

Sterilisation & Packaging

Supplied clean and non-sterile for hospital autoclaving, or gamma / ETO sterilised on request. Individually barcoded and lot-traceable with a validated 5-year shelf life. Matching instrument sets and loan kits available for all listed sizes.

Available for export to 40+ countries. Contact our team for full catalogue, surgical technique guide and territory pricing.

Certifications

- ISO 13485
- CE Mark

Contact

- **Address:** 14 Harley Medical Court, London W1G 9PF, United Kingdom
- **Telephone:** [+91 9319821922](tel:+919319821922)
- **Email:** exports@orthomedimplants.com

This document was generated from <https://www.jksigma.in/implants/calcaneal-locking-plate/>