

BONE SCREWS

Cancellous Bone Screw

Product code: SCR-CANC



Specifications

Specification	Value
Material	Titanium Ti-6Al-4V
Diameter range	4.0mm / 6.5mm
Length range	20mm – 120mm
Thread type	Coarse cancellous, partial/full
Availability	In stock — ready to ship

Available sizes and configurations

- 4.0mm
- 6.5mm

Product description

Coarse-thread cancellous screws with deep threads for maximum purchase in metaphyseal and epiphyseal cancellous bone.

Key Features

- Deep coarse thread
- 16mm and 32mm thread lengths
- Partial thread for lag effect
- Washer compatible

Technical Specifications

Product Code	SCR-CANC
Material	Titanium Ti-6Al-4V
Diameter	4.0mm / 6.5mm
Length Range	20mm – 120mm
Available Sizes	4.0mm, 6.5mm
Thread / Locking	Coarse cancellous, partial/full
Certifications	ISO 13485, CE Mark

Sterilisation & Packaging

Supplied clean and non-sterile for hospital autoclaving, or gamma / ETO sterilised on request. Individually barcoded and lot-traceable with a validated 5-year shelf life. Matching instrument sets and loan kits available for all listed sizes.

Available for export to 40+ countries. Contact our team for full catalogue, surgical technique guide and territory pricing.

Certifications

- ISO 13485
- CE Mark

Contact

- **Address:** 14 Harley Medical Court, London W1G 9PF, United Kingdom
- **Telephone:** [+91 9319821922](tel:+919319821922)
- **Email:** exports@orthomedimplants.com

This document was generated from <https://www.jksigma.in/implants/cancellous-bone-screw/>