

BONE SCREWS

Cannulated Cancellous Screw

Product code: SCR-CANN



Specifications

Specification	Value
Material	Titanium Ti-6Al-4V
Diameter range	3.5mm – 7.3mm
Length range	20mm – 130mm
Thread type	Cannulated, partial thread
Availability	In stock — ready to ship

Available sizes and configurations

- 3.5mm
- 4.5mm
- 6.5mm
- 7.3mm

Product description

Guidewire-cannulated cancellous screws for percutaneous, accurately guided fixation of femoral neck and periarticular fractures.

Key Features

- Guidewire-assisted placement
- Partial thread creates interfragmentary compression
- Self-drilling / self-tapping tip
- Washers available

Technical Specifications

Product Code	SCR-CANN
Material	Titanium Ti-6Al-4V
Diameter	3.5mm – 7.3mm
Length Range	20mm – 130mm
Available Sizes	3.5mm, 4.5mm, 6.5mm, 7.3mm
Thread / Locking	Cannulated, partial thread
Certifications	ISO 13485, CE Mark, FDA 510(k)

Sterilisation & Packaging

Supplied clean and non-sterile for hospital autoclaving, or gamma / ETO sterilised on request. Individually barcoded and lot-traceable with a validated 5-year shelf life. Matching instrument sets and loan kits available for all listed sizes.

Available for export to 40+ countries. Contact our team for full catalogue, surgical technique guide and territory pricing.

Certifications

- ISO 13485
- CE Mark
- FDA 510(k)

Contact

- **Address:** 14 Harley Medical Court, London W1G 9PF, United Kingdom
- **Telephone:** [+91 9319821922](tel:+919319821922)
- **Email:** exports@orthomedimplants.com

This document was generated from <https://www.jksigma.in/implants/cannulated-cancellous-screw/>