


PLATES & FIXATORS

Clavicle Locking Plate

Product code: CLAV-PL



Specifications

Specification	Value
Material	Titanium Ti-6Al-4V
Diameter range	2.7mm / 3.5mm screws
Length range	80mm – 140mm
Thread type	Locking, pre-contoured
Availability	In stock — ready to ship

Available sizes and configurations

- 6
- 7
- 8 holes

Product description

A superior, anatomically pre-contoured clavicle plate for displaced midshaft and lateral clavicle fractures, reducing the need for intra-operative bending.

Key Features

- S-shaped anatomic pre-contour
- Left/right specific
- Lateral extension for distal fractures
- Low-profile locking heads

Technical Specifications

Product Code	CLAV-PL
Material	Titanium Ti-6Al-4V
Diameter	2.7mm / 3.5mm screws
Length Range	80mm – 140mm
Available Sizes	6, 7, 8 holes
Thread / Locking	Locking, pre-contoured
Certifications	ISO 13485, CE Mark

Sterilisation & Packaging

Supplied clean and non-sterile for hospital autoclaving, or gamma / ETO sterilised on request. Individually barcoded and lot-traceable with a validated 5-year shelf life. Matching instrument sets and loan kits available for all listed sizes.

Available for export to 40+ countries. Contact our team for full catalogue, surgical technique guide and territory pricing.

Certifications

- ISO 13485
- CE Mark

Contact

- **Address:** 14 Harley Medical Court, London W1G 9PF, United Kingdom
- **Telephone:** [+91 9319821922](tel:+919319821922)
- **Email:** exports@orthomedimplants.com

This document was generated from <https://www.jksigma.in/implants/clavicle-locking-plate/>