


BONE SCREWS

Cortical Bone Screw

Product code: SCR-CORT



Specifications

Specification	Value
Material	Titanium Ti-6Al-4V
Diameter range	2.7mm – 4.5mm
Length range	10mm – 110mm
Thread type	Fine cortical, self-tapping
Availability	In stock — ready to ship

Available sizes and configurations

- 2.7mm
- 3.5mm
- 4.5mm

Product description

Fully threaded self-tapping cortical screws for diaphyseal fixation and plate attachment in dense cortical bone.

Key Features

- Fine-pitch cortical thread
- Self-tapping tip
- Hex / star recess drive
- Full length range in 2mm steps

Technical Specifications

Product Code	SCR-CORT
Material	Titanium Ti-6Al-4V
Diameter	2.7mm – 4.5mm
Length Range	10mm – 110mm
Available Sizes	2.7mm, 3.5mm, 4.5mm
Thread / Locking	Fine cortical, self-tapping
Certifications	ISO 13485, CE Mark, FDA 510(k)

Sterilisation & Packaging

Supplied clean and non-sterile for hospital autoclaving, or gamma / ETO sterilised on request. Individually barcoded and lot-traceable with a validated 5-year shelf life. Matching instrument sets and loan kits available for all listed sizes.

Available for export to 40+ countries. Contact our team for full catalogue, surgical technique guide and territory pricing.

Certifications

- ISO 13485
- CE Mark
- FDA 510(k)

Contact

- **Address:** 14 Harley Medical Court, London W1G 9PF, United Kingdom
- **Telephone:** [+91 9319821922](tel:+919319821922)
- **Email:** exports@orthomedimplants.com

This document was generated from <https://www.jksigma.in/implants/cortical-bone-screw/>