


PLATES & FIXATORS

Distal Femur Locking Plate (LISS)

Product code: DFP-LISS



Specifications

Specification	Value
Material	Titanium Ti-6Al-4V
Diameter range	5.0mm locking screws
Length range	160mm – 360mm
Thread type	Fixed-angle locking
Availability	In stock — ready to ship

Available sizes and configurations

- 5
- 9
- 13 holes

Product description

A lateral distal femur plate designed for minimally invasive submuscular insertion (MIPO), acting as an internal fixator for comminuted supracondylar fractures.

Key Features

- Anatomic lateral contour, left/right
- Percutaneous aiming block
- Fixed-angle condylar locking screws
- Preserves periosteal blood supply

Technical Specifications

Product Code	DFP-LISS
Material	Titanium Ti-6Al-4V
Diameter	5.0mm locking screws
Length Range	160mm – 360mm
Available Sizes	5, 9, 13 holes
Thread / Locking	Fixed-angle locking
Certifications	ISO 13485, CE Mark

Sterilisation & Packaging

Supplied clean and non-sterile for hospital autoclaving, or gamma / ETO sterilised on request. Individually barcoded and lot-traceable with a validated 5-year shelf life. Matching instrument sets and loan kits available for all listed sizes.

Available for export to 40+ countries. Contact our team for full catalogue, surgical technique guide and territory pricing.

Certifications

- ISO 13485
- CE Mark

Contact

- **Address:** 14 Harley Medical Court, London W1G 9PF, United Kingdom
- **Telephone:** [+91 9319821922](tel:+919319821922)
- **Email:** exports@orthomedimplants.com

This document was generated from <https://www.jksigma.in/implants/distal-femur-locking-plate-liss/>