


PLATES & FIXATORS

Dynamic Compression Plate (DCP)

Product code: DCP



Specifications

Specification	Value
Material	Stainless Steel 316L
Diameter range	3.5mm / 4.5mm screws
Length range	55mm – 290mm
Thread type	Inclined-cylinder holes
Availability	In stock — ready to ship

Available sizes and configurations

- 4–16 holes

Product description

A classic dynamic compression plate whose inclined oval holes generate axial compression as

eccentrically placed screws are tightened.

Key Features

- Spherical-gliding hole principle
- Up to 1mm compression per hole
- Narrow and broad variants
- Proven AO fixation standard

Technical Specifications

Product Code	DCP
Material	Stainless Steel 316L
Diameter	3.5mm / 4.5mm screws
Length Range	55mm – 290mm
Available Sizes	4–16 holes
Thread / Locking	Inclined-cylinder holes
Certifications	ISO 13485, CE Mark

Sterilisation & Packaging

Supplied clean and non-sterile for hospital autoclaving, or gamma / ETO sterilised on request. Individually barcoded and lot-traceable with a validated 5-year shelf life. Matching instrument sets and loan kits available for all listed sizes.

Available for export to 40+ countries. Contact our team for full catalogue, surgical technique guide and territory pricing.

Certifications

- ISO 13485
- CE Mark

Contact

- **Address:** 14 Harley Medical Court, London W1G 9PF, United Kingdom

- **Telephone:** [+91 9319821922](tel:+919319821922)
- **Email:** exports@orthomedimplants.com

This document was generated from <https://www.jksigma.in/implants/dynamic-compression-plate-dcp/>