


PLATES & FIXATORS

Dynamic Condylar Screw (DCS) Plate

Product code: DCS



Specifications

Specification	Value
Material	Stainless Steel 316L
Diameter range	95° angle
Length range	90mm – 250mm
Thread type	Lag screw + side plate
Availability	Available on order

Available sizes and configurations

- 4–12 holes

Product description

A fixed 95° condylar screw-and-plate construct for distal femoral and subtrochanteric fractures requiring

rigid angular stability.

Key Features

- Fixed 95° blade-plate geometry
- Large condylar lag screw
- Handles supracondylar comminution
- Compression via articulated tensioner

Technical Specifications

Product Code	DCS
Material	Stainless Steel 316L
Diameter	95° angle
Length Range	90mm – 250mm
Available Sizes	4–12 holes
Thread / Locking	Lag screw + side plate
Certifications	ISO 13485, CE Mark

Sterilisation & Packaging

Supplied clean and non-sterile for hospital autoclaving, or gamma / ETO sterilised on request. Individually barcoded and lot-traceable with a validated 5-year shelf life. Matching instrument sets and loan kits available for all listed sizes.

Available for export to 40+ countries. Contact our team for full catalogue, surgical technique guide and territory pricing.

Certifications

- ISO 13485
- CE Mark

Contact

- **Address:** 14 Harley Medical Court, London W1G 9PF, United Kingdom

- **Telephone:** [+91 9319821922](tel:+919319821922)
- **Email:** exports@orthomedimplants.com

This document was generated from <https://www.jksigma.in/implants/dynamic-condylar-screw-dcs-plate/>