


PLATES & FIXATORS

Dynamic Hip Screw (DHS) Plate

Product code: DHS



Specifications

Specification	Value
Material	Stainless Steel 316L
Diameter range	135°, 140°, 145° barrel
Length range	38mm – 200mm side plate
Thread type	Lag screw + barrel plate
Availability	In stock — ready to ship

Available sizes and configurations

- 2–8 holes

Product description

The gold-standard sliding hip screw for stable intertrochanteric and femoral neck fractures, allowing

controlled fracture impaction under load.

Key Features

- 135° / 140° / 145° barrel angles
- Controlled lag-screw sliding
- 2 to 8 hole side plates
- Optional anti-rotation trochanteric screw

Technical Specifications

Product Code	DHS
Material	Stainless Steel 316L
Diameter	135°, 140°, 145° barrel
Length Range	38mm – 200mm side plate
Available Sizes	2–8 holes
Thread / Locking	Lag screw + barrel plate
Certifications	ISO 13485, CE Mark, FDA 510(k)

Sterilisation & Packaging

Supplied clean and non-sterile for hospital autoclaving, or gamma / ETO sterilised on request. Individually barcoded and lot-traceable with a validated 5-year shelf life. Matching instrument sets and loan kits available for all listed sizes.

Available for export to 40+ countries. Contact our team for full catalogue, surgical technique guide and territory pricing.

Certifications

- ISO 13485
- CE Mark
- FDA 510(k)

Contact

- **Address:** 14 Harley Medical Court, London W1G 9PF, United Kingdom
- **Telephone:** [+91 9319821922](tel:+919319821922)
- **Email:** exports@orthomedimplants.com

This document was generated from <https://www.jksigma.in/implants/dynamic-hip-screw-dhs-plate/>