

INTRAMEDULLARY NAILS

Expert Suprapatellar Tibial Nail

Product code: IMN-TIB-SP



Specifications

Specification	Value
Material	Titanium Ti-6Al-4V
Diameter range	8.0mm – 11.0mm
Length range	255mm – 420mm
Thread type	Cannulated
Availability	In stock — ready to ship

Available sizes and configurations

- 8mm
- 9mm
- 10mm
- 11mm

Product description

Optimised for semi-extended suprapatellar insertion, extending intramedullary fixation to very proximal and very distal tibial fractures.

Key Features

- Five proximal + four distal locking holes
- Angular stable locking bushings
- Semi-extended positioning reduces malalignment
- Radiolucent instrumentation

Technical Specifications

Product Code	IMN-TIB-SP
Material	Titanium Ti-6Al-4V
Diameter	8.0mm – 11.0mm
Length Range	255mm – 420mm
Available Sizes	8mm, 9mm, 10mm, 11mm
Thread / Locking	Cannulated
Certifications	ISO 13485, CE Mark

Sterilisation & Packaging

Supplied clean and non-sterile for hospital autoclaving, or gamma / ETO sterilised on request. Individually barcoded and lot-traceable with a validated 5-year shelf life. Matching instrument sets and loan kits available for all listed sizes.

Available for export to 40+ countries. Contact our team for full catalogue, surgical technique guide and territory pricing.

Certifications

- ISO 13485
- CE Mark

Contact

- **Address:** 14 Harley Medical Court, London W1G 9PF, United Kingdom
- **Telephone:** [+91 9319821922](tel:+919319821922)
- **Email:** exports@orthomedimplants.com

This document was generated from <https://www.jksigma.in/implants/expert-suprapatellar-tibial-nail/>