


INTRAMEDULLARY NAILS

Gamma-Type Trochanteric Nail

Product code: IMN-TROCH



Specifications

Specification	Value
Material	Stainless Steel 316L
Diameter range	10.0mm – 14.0mm
Length range	180mm – 420mm
Thread type	Lag screw
Availability	Available on order

Available sizes and configurations

- 10mm
- 11mm
- 12mm
- 14mm

Product description

A short trochanteric nail with a large-diameter lag screw for stable and unstable intertrochanteric fractures, allowing early weight bearing.

Key Features

- Large 10.5mm lag screw
- Set screw controls sliding/rotation
- Short and long lengths
- Single-shot distal locking jig

Technical Specifications

Product Code	IMN-TROCH
Material	Stainless Steel 316L
Diameter	10.0mm – 14.0mm
Length Range	180mm – 420mm
Available Sizes	10mm, 11mm, 12mm, 14mm
Thread / Locking	Lag screw
Certifications	ISO 13485, CE Mark

Sterilisation & Packaging

Supplied clean and non-sterile for hospital autoclaving, or gamma / ETO sterilised on request. Individually barcoded and lot-traceable with a validated 5-year shelf life. Matching instrument sets and loan kits available for all listed sizes.

Available for export to 40+ countries. Contact our team for full catalogue, surgical technique guide and territory pricing.

Certifications

- ISO 13485
- CE Mark

Contact

- **Address:** 14 Harley Medical Court, London W1G 9PF, United Kingdom
- **Telephone:** [+91 9319821922](tel:+919319821922)
- **Email:** exports@orthomedimplants.com

This document was generated from <https://www.jksigma.in/implants/gamma-type-trochanteric-nail/>