


BONE SCREWS

Headless Compression Screw

Product code: SCR-HCS



Specifications

Specification	Value
Material	Titanium Ti-6Al-4V
Diameter range	2.5mm – 4.5mm
Length range	8mm – 60mm
Thread type	Dual-pitch, fully threaded
Availability	In stock — ready to ship

Available sizes and configurations

- 2.5mm
- 3.0mm
- 4.5mm

Product description

A dual-pitch headless screw that buries beneath the articular surface and generates continuous compression, ideal for scaphoid and small-fragment fractures.

Key Features

- Variable dual-pitch thread
- Fully buried, no prominent head
- Continuous interfragmentary compression
- Scaphoid, radial head, OCD fixation

Technical Specifications

Product Code	SCR-HCS
Material	Titanium Ti-6Al-4V
Diameter	2.5mm – 4.5mm
Length Range	8mm – 60mm
Available Sizes	2.5mm, 3.0mm, 4.5mm
Thread / Locking	Dual-pitch, fully threaded
Certifications	ISO 13485, CE Mark

Sterilisation & Packaging

Supplied clean and non-sterile for hospital autoclaving, or gamma / ETO sterilised on request. Individually barcoded and lot-traceable with a validated 5-year shelf life. Matching instrument sets and loan kits available for all listed sizes.

Available for export to 40+ countries. Contact our team for full catalogue, surgical technique guide and territory pricing.

Certifications

- ISO 13485
- CE Mark

Contact

- **Address:** 14 Harley Medical Court, London W1G 9PF, United Kingdom
- **Telephone:** [+91 9319821922](tel:+919319821922)
- **Email:** exports@orthomedimplants.com

This document was generated from <https://www.jksigma.in/implants/headless-compression-screw/>