


INTRAMEDULLARY NAILS

Humerus Interlocking Nail

Product code: IMN-HUM



Specifications

Specification	Value
Material	Titanium Ti-6Al-4V
Diameter range	6.0mm – 9.0mm
Length range	160mm – 300mm
Thread type	Cannulated
Availability	In stock — ready to ship

Available sizes and configurations

- 6mm
- 7mm
- 8mm
- 9mm

Product description

A low-profile humeral nail for proximal, shaft and combined humeral fractures, with a rounded proximal end that minimises rotator-cuff irritation.

Key Features

- Antegrade and retrograde insertion
- Multi-directional proximal locking
- Cuff-friendly rounded proximal tip
- Suture holes for tuberosity fixation

Technical Specifications

Product Code	IMN-HUM
Material	Titanium Ti-6Al-4V
Diameter	6.0mm – 9.0mm
Length Range	160mm – 300mm
Available Sizes	6mm, 7mm, 8mm, 9mm
Thread / Locking	Cannulated
Certifications	ISO 13485, CE Mark

Sterilisation & Packaging

Supplied clean and non-sterile for hospital autoclaving, or gamma / ETO sterilised on request. Individually barcoded and lot-traceable with a validated 5-year shelf life. Matching instrument sets and loan kits available for all listed sizes.

Available for export to 40+ countries. Contact our team for full catalogue, surgical technique guide and territory pricing.

Certifications

- ISO 13485
- CE Mark

Contact

- **Address:** 14 Harley Medical Court, London W1G 9PF, United Kingdom
- **Telephone:** [+91 9319821922](tel:+919319821922)
- **Email:** exports@orthomedimplants.com

This document was generated from <https://www.jksigma.in/implants/humerus-interlocking-nail/>