


PLATES & FIXATORS

Locking Compression Plate — Broad

Product code: LCP-BRD



Specifications

Specification	Value
Material	Titanium Ti-6Al-4V
Diameter range	4.5mm / 5.0mm screws
Length range	110mm – 340mm
Thread type	Combi-hole (staggered)
Availability	In stock — ready to ship

Available sizes and configurations

- 6–18 holes

Product description

A broad LCP with staggered combi-holes for large-bone diaphyseal fractures requiring maximum

construct strength, such as the femur and humerus.

Key Features

- Staggered holes reduce stress risers
- High bending stiffness for femoral shaft
- Locking + compression in one plate
- 6 to 18 hole lengths

Technical Specifications

Product Code	LCP-BRD
Material	Titanium Ti-6Al-4V
Diameter	4.5mm / 5.0mm screws
Length Range	110mm – 340mm
Available Sizes	6–18 holes
Thread / Locking	Combi-hole (staggered)
Certifications	ISO 13485, CE Mark, FDA 510(k)

Sterilisation & Packaging

Supplied clean and non-sterile for hospital autoclaving, or gamma / ETO sterilised on request. Individually barcoded and lot-traceable with a validated 5-year shelf life. Matching instrument sets and loan kits available for all listed sizes.

Available for export to 40+ countries. Contact our team for full catalogue, surgical technique guide and territory pricing.

Certifications

- ISO 13485
- CE Mark
- FDA 510(k)

Contact

- **Address:** 14 Harley Medical Court, London W1G 9PF, United Kingdom
- **Telephone:** [+91 9319821922](tel:+919319821922)
- **Email:** exports@orthomedimplants.com

This document was generated from <https://www.jksigma.in/implants/locking-compression-plate-broad/>