


PLATES & FIXATORS

Locking Compression Plate — Narrow

Product code: LCP-NAR



Specifications

Specification	Value
Material	Titanium Ti-6Al-4V
Diameter range	3.5mm / 4.5mm screws
Length range	56mm – 250mm
Thread type	Combi-hole (locking + compression)
Availability	In stock — ready to ship

Available sizes and configurations

- 4–14 holes

Product description

The workhorse locking compression plate with combi-holes that accept both locking and standard

cortical screws in the same construct for diaphyseal fractures.

Key Features

- Combi-holes: locking OR compression per hole
- Low-profile, limited-contact undersurface
- Colour-anodised for size identification
- 4 to 14 hole lengths

Technical Specifications

Product Code	LCP-NAR
Material	Titanium Ti-6Al-4V
Diameter	3.5mm / 4.5mm screws
Length Range	56mm – 250mm
Available Sizes	4–14 holes
Thread / Locking	Combi-hole (locking + compression)
Certifications	ISO 13485, CE Mark, FDA 510(k)

Sterilisation & Packaging

Supplied clean and non-sterile for hospital autoclaving, or gamma / ETO sterilised on request. Individually barcoded and lot-traceable with a validated 5-year shelf life. Matching instrument sets and loan kits available for all listed sizes.

Available for export to 40+ countries. Contact our team for full catalogue, surgical technique guide and territory pricing.

Certifications

- ISO 13485
- CE Mark
- FDA 510(k)

Contact

- **Address:** 14 Harley Medical Court, London W1G 9PF, United Kingdom
- **Telephone:** [+91 9319821922](tel:+919319821922)
- **Email:** exports@orthomedimplants.com

This document was generated from <https://www.jksigma.in/implants/locking-compression-plate-narrow/>