


PLATES & FIXATORS

Monolateral External Fixator

Product code: EXFIX-MONO



Specifications

Specification	Value
Material	Aluminium + Stainless Steel
Diameter range	4mm / 5mm / 6mm pins
Length range	100mm – 350mm body
Thread type	Schanz pins
Availability	In stock — ready to ship

Available sizes and configurations

- Small
- Standard
- Large

Product description

A rigid modular monolateral frame for temporary damage-control fixation and definitive treatment of open tibial and femoral fractures.

Key Features

- Radiolucent carbon-fibre body option
- Multi-planar pin clamps
- Compression/distraction unit
- Reusable, autoclavable components

Technical Specifications

Product Code	EXFIX-MONO
Material	Aluminium + Stainless Steel
Diameter	4mm / 5mm / 6mm pins
Length Range	100mm – 350mm body
Available Sizes	Small, Standard, Large
Thread / Locking	Schanz pins
Certifications	ISO 13485, CE Mark

Sterilisation & Packaging

Supplied clean and non-sterile for hospital autoclaving, or gamma / ETO sterilised on request. Individually barcoded and lot-traceable with a validated 5-year shelf life. Matching instrument sets and loan kits available for all listed sizes.

Available for export to 40+ countries. Contact our team for full catalogue, surgical technique guide and territory pricing.

Certifications

- ISO 13485
- CE Mark

Contact

- **Address:** 14 Harley Medical Court, London W1G 9PF, United Kingdom
- **Telephone:** [+91 9319821922](tel:+919319821922)
- **Email:** exports@orthomedimplants.com

This document was generated from <https://www.jksigma.in/implants/monolateral-external-fixator/>