


PLATES & FIXATORS

One-Third Tubular Plate

Product code: TUB-13



Specifications

Specification	Value
Material	Stainless Steel 316L
Diameter range	3.5mm screws
Length range	35mm – 150mm
Thread type	Oval compression holes
Availability	In stock — ready to ship

Available sizes and configurations

- 3–12 holes

Product description

A thin, low-profile plate for lateral malleolus, olecranon and other subcutaneous locations where

soft-tissue coverage is limited.

Key Features

- 1.0mm thickness, low profile
- Collar hole for lag-screw compression
- Easily contoured
- Ideal for lateral malleolus

Technical Specifications

Product Code	TUB-13
Material	Stainless Steel 316L
Diameter	3.5mm screws
Length Range	35mm – 150mm
Available Sizes	3–12 holes
Thread / Locking	Oval compression holes
Certifications	ISO 13485, CE Mark

Sterilisation & Packaging

Supplied clean and non-sterile for hospital autoclaving, or gamma / ETO sterilised on request. Individually barcoded and lot-traceable with a validated 5-year shelf life. Matching instrument sets and loan kits available for all listed sizes.

Available for export to 40+ countries. Contact our team for full catalogue, surgical technique guide and territory pricing.

Certifications

- ISO 13485
- CE Mark

Contact

- **Address:** 14 Harley Medical Court, London W1G 9PF, United Kingdom

- **Telephone:** [+91 9319821922](tel:+919319821922)
- **Email:** exports@orthomedimplants.com

This document was generated from <https://www.jksigma.in/implants/one-third-tubular-plate/>