


INTRAMEDULLARY NAILS

Proximal Femoral Nail Anti-Rotation (PFNA)

Product code: IMN-PFNA



Specifications

Specification	Value
Material	Titanium Ti-6Al-4V
Diameter range	10.0mm – 12.0mm
Length range	240mm – 420mm
Thread type	Helical blade
Availability	In stock — ready to ship

Available sizes and configurations

- 10mm
- 11mm
- 12mm

Product description

A single helical-blade cephalomedullary nail that compacts cancellous bone for superior purchase in osteoporotic pertrochanteric fractures.

Key Features

- Helical blade compacts, not cuts, bone
- Rotational and angular stability in one step
- Ideal for osteoporotic geriatric fractures
- 6° mediolateral angle for lateral entry

Technical Specifications

Product Code	IMN-PFNA
Material	Titanium Ti-6Al-4V
Diameter	10.0mm – 12.0mm
Length Range	240mm – 420mm
Available Sizes	10mm, 11mm, 12mm
Thread / Locking	Helical blade
Certifications	ISO 13485, CE Mark, FDA 510(k)

Sterilisation & Packaging

Supplied clean and non-sterile for hospital autoclaving, or gamma / ETO sterilised on request. Individually barcoded and lot-traceable with a validated 5-year shelf life. Matching instrument sets and loan kits available for all listed sizes.

Available for export to 40+ countries. Contact our team for full catalogue, surgical technique guide and territory pricing.

Certifications

- ISO 13485
- CE Mark
- FDA 510(k)

Contact

- **Address:** 14 Harley Medical Court, London W1G 9PF, United Kingdom
- **Telephone:** [+91 9319821922](tel:+919319821922)
- **Email:** exports@orthomedimplants.com

This document was generated from <https://www.jksigma.in/implants/proximal-femoral-nail-antirotation-pfna/>