


INTRAMEDULLARY NAILS

Proximal Femoral Nail (PFN)

Product code: IMN-PFN



Specifications

Specification	Value
Material	Titanium Ti-6Al-4V
Diameter range	9.0mm – 12.0mm
Length range	240mm – 420mm
Thread type	Cannulated
Availability	In stock — ready to ship

Available sizes and configurations

- 9mm
- 10mm
- 11mm
- 12mm

Product description

A two-screw cephalomedullary nail for pertrochanteric and subtrochanteric femoral fractures, providing controlled compression across the fracture with anti-rotation support.

Key Features

- 130° / 135° neck-shaft angle options
- Neck screw plus anti-rotation pin
- Distal static and dynamic locking
- Valgus-reducing proximal geometry

Technical Specifications

Product Code	IMN-PFN
Material	Titanium Ti-6Al-4V
Diameter	9.0mm – 12.0mm
Length Range	240mm – 420mm
Available Sizes	9mm, 10mm, 11mm, 12mm
Thread / Locking	Cannulated
Certifications	ISO 13485, CE Mark, FDA 510(k)

Sterilisation & Packaging

Supplied clean and non-sterile for hospital autoclaving, or gamma / ETO sterilised on request. Individually barcoded and lot-traceable with a validated 5-year shelf life. Matching instrument sets and loan kits available for all listed sizes.

Available for export to 40+ countries. Contact our team for full catalogue, surgical technique guide and territory pricing.

Certifications

- ISO 13485
- CE Mark
- FDA 510(k)

Contact

- **Address:** 14 Harley Medical Court, London W1G 9PF, United Kingdom
- **Telephone:** [+91 9319821922](tel:+919319821922)
- **Email:** exports@orthomedimplants.com

This document was generated from <https://www.jksigma.in/implants/proximal-femoral-nail-pfn/>