


PLATES & FIXATORS

Proximal Humerus Locking Plate (PHILOS)

Product code: PHP-PHILOS



Specifications

Specification	Value
Material	Titanium Ti-6Al-4V
Diameter range	3.5mm locking screws
Length range	90mm – 220mm
Thread type	Multidirectional locking
Availability	In stock — ready to ship

Available sizes and configurations

- 3
- 5
- 8 holes

Product description

An anatomically contoured proximal humerus plate with a divergent locking-screw array and suture holes for two-, three- and four-part fractures.

Key Features

- Up to 9 divergent head locking screws
- Calcar screw supports medial column
- Suture holes for rotator-cuff fixation
- Low-profile anatomic contour

Technical Specifications

Product Code	PHP-PHILOS
Material	Titanium Ti-6Al-4V
Diameter	3.5mm locking screws
Length Range	90mm – 220mm
Available Sizes	3, 5, 8 holes
Thread / Locking	Multidirectional locking
Certifications	ISO 13485, CE Mark, FDA 510(k)

Sterilisation & Packaging

Supplied clean and non-sterile for hospital autoclaving, or gamma / ETO sterilised on request. Individually barcoded and lot-traceable with a validated 5-year shelf life. Matching instrument sets and loan kits available for all listed sizes.

Available for export to 40+ countries. Contact our team for full catalogue, surgical technique guide and territory pricing.

Certifications

- ISO 13485
- CE Mark
- FDA 510(k)

Contact

- **Address:** 14 Harley Medical Court, London W1G 9PF, United Kingdom
- **Telephone:** [+91 9319821922](tel:+919319821922)
- **Email:** exports@orthomedimplants.com

This document was generated from <https://www.jksigma.in/implants/proximal-humerus-locking-plate-philos/>