

PLATES & FIXATORS

Proximal Tibia Locking Plate

Product code: PTP



Specifications

Specification	Value
Material	Titanium Ti-6Al-4V
Diameter range	3.5mm / 4.5mm screws
Length range	110mm – 290mm
Thread type	Locking + compression
Availability	In stock — ready to ship

Available sizes and configurations

- 4
- 6
- 10 holes

Product description

An anatomic lateral tibial plateau plate supporting the articular surface with a raft of subchondral locking screws for split-depression fractures.

Key Features

- Subchondral raft screw configuration
- Anatomic lateral plateau contour
- Left/right specific
- Combi-holes in the shaft

Technical Specifications

Product Code	PTP
Material	Titanium Ti-6Al-4V
Diameter	3.5mm / 4.5mm screws
Length Range	110mm – 290mm
Available Sizes	4, 6, 10 holes
Thread / Locking	Locking + compression
Certifications	ISO 13485, CE Mark

Sterilisation & Packaging

Supplied clean and non-sterile for hospital autoclaving, or gamma / ETO sterilised on request. Individually barcoded and lot-traceable with a validated 5-year shelf life. Matching instrument sets and loan kits available for all listed sizes.

Available for export to 40+ countries. Contact our team for full catalogue, surgical technique guide and territory pricing.

Certifications

- ISO 13485
- CE Mark

Contact

- **Address:** 14 Harley Medical Court, London W1G 9PF, United Kingdom
- **Telephone:** [+91 9319821922](tel:+919319821922)
- **Email:** exports@orthomedimplants.com

This document was generated from <https://www.jksigma.in/implants/proximal-tibia-locking-plate/>