


INTRAMEDULLARY NAILS

# Reconstruction Femoral Nail

Product code: IMN-RECON



## Specifications

| Specification  | Value                    |
|----------------|--------------------------|
| Material       | Titanium Ti-6Al-4V       |
| Diameter range | 10.0mm – 13.0mm          |
| Length range   | 300mm – 460mm            |
| Thread type    | Dual proximal screws     |
| Availability   | In stock — ready to ship |

## Available sizes and configurations

- 10mm
- 11mm
- 12mm
- 13mm

## Product description

A reconstruction nail with two cephalic screws directed into the femoral neck for ipsilateral neck-and-shaft fractures and pathological lesions.

### Key Features

- Two divergent cephalic locking screws
- Handles ipsilateral neck + shaft fractures
- Prophylactic fixation of metastatic lesions
- Cannulated full length

### Technical Specifications

| Product Code     | IMN-RECON              |
|------------------|------------------------|
| Material         | Titanium Ti-6Al-4V     |
| Diameter         | 10.0mm – 13.0mm        |
| Length Range     | 300mm – 460mm          |
| Available Sizes  | 10mm, 11mm, 12mm, 13mm |
| Thread / Locking | Dual proximal screws   |
| Certifications   | ISO 13485, CE Mark     |

### Sterilisation & Packaging

Supplied clean and non-sterile for hospital autoclaving, or gamma / ETO sterilised on request. Individually barcoded and lot-traceable with a validated 5-year shelf life. Matching instrument sets and loan kits available for all listed sizes.

*Available for export to 40+ countries. Contact our team for full catalogue, surgical technique guide and territory pricing.*

### Certifications

- ISO 13485
- CE Mark

## Contact

- **Address:** 14 Harley Medical Court, London W1G 9PF, United Kingdom
- **Telephone:** [+91 9319821922](tel:+919319821922)
- **Email:** [exports@orthomedimplants.com](mailto:exports@orthomedimplants.com)

---

This document was generated from <https://www.jksigma.in/implants/reconstruction-femoral-nail/>