


PLATES & FIXATORS

Reconstruction Plate

Product code: RECON-PL



Specifications

Specification	Value
Material	Stainless Steel 316L
Diameter range	3.5mm screws
Length range	70mm – 300mm
Thread type	Notched, contourable
Availability	In stock — ready to ship

Available sizes and configurations

- 6–20 holes

Product description

A notched plate that can be bent and twisted in three planes for pelvic, acetabular, clavicle and

periarticular fractures with complex geometry.

Key Features

- Contourable in three planes
- Deep notches for tight bending radii
- Straight and pre-curved variants
- 3.5mm cortical / locking compatible

Technical Specifications

Product Code	RECON-PL
Material	Stainless Steel 316L
Diameter	3.5mm screws
Length Range	70mm – 300mm
Available Sizes	6–20 holes
Thread / Locking	Notched, contourable
Certifications	ISO 13485, CE Mark

Sterilisation & Packaging

Supplied clean and non-sterile for hospital autoclaving, or gamma / ETO sterilised on request. Individually barcoded and lot-traceable with a validated 5-year shelf life. Matching instrument sets and loan kits available for all listed sizes.

Available for export to 40+ countries. Contact our team for full catalogue, surgical technique guide and territory pricing.

Certifications

- ISO 13485
- CE Mark

Contact

- **Address:** 14 Harley Medical Court, London W1G 9PF, United Kingdom

- **Telephone:** [+91 9319821922](tel:+919319821922)
- **Email:** exports@orthomedimplants.com

This document was generated from <https://www.jksigma.in/implants/reconstruction-plate/>