


INTRAMEDULLARY NAILS

Retrograde Femoral Nail

Product code: IMN-FEM-RG



Specifications

Specification	Value
Material	Titanium Ti-6Al-4V
Diameter range	10.0mm – 13.0mm
Length range	160mm – 400mm
Thread type	Cannulated
Availability	In stock — ready to ship

Available sizes and configurations

- 10mm
- 11mm
- 12mm
- 13mm

Product description

Designed for retrograde insertion through the intercondylar notch, ideal for distal femoral, periprosthetic and polytrauma fractures where antegrade access is unsuitable.

Key Features

- Short and long variants on one platform
- Spiral blade and screw distal locking
- Recessed end cap prevents knee impingement
- Oblique and AP distal locking holes

Technical Specifications

Product Code	IMN-FEM-RG
Material	Titanium Ti-6Al-4V
Diameter	10.0mm – 13.0mm
Length Range	160mm – 400mm
Available Sizes	10mm, 11mm, 12mm, 13mm
Thread / Locking	Cannulated
Certifications	ISO 13485, CE Mark

Sterilisation & Packaging

Supplied clean and non-sterile for hospital autoclaving, or gamma / ETO sterilised on request. Individually barcoded and lot-traceable with a validated 5-year shelf life. Matching instrument sets and loan kits available for all listed sizes.

Available for export to 40+ countries. Contact our team for full catalogue, surgical technique guide and territory pricing.

Certifications

- ISO 13485
- CE Mark

Contact

- **Address:** 14 Harley Medical Court, London W1G 9PF, United Kingdom
- **Telephone:** [+91 9319821922](tel:+919319821922)
- **Email:** exports@orthomedimplants.com

This document was generated from <https://www.jksigma.in/implants/retrograde-femoral-nail/>