

Manufacturing

World-Class Manufacturing for Life-Critical Implants

Raw Material

Ti-6Al-4V ELI and 316L SS sourced exclusively from certified aerospace/medical mills. Each billet arrives with a Mill Test Certificate; our incoming inspection team performs independent hardness and spectroscopic verification before any material enters machining.

CNC Machining

28 five-axis CNC centres with on-machine probing for real-time dimensional feedback. Toolpaths generated from validated CAM programmes under document control. Critical features held to ≤ 0.01 mm tolerance.

Surface Treatment

- **Electropolishing** — $R_a \leq 0.05$ μ m on articulating surfaces
- **Type III Hard Anodising** — HV 600–800 for titanium implants
- **HA Plasma Spray** — 50–75 μ m coating, $\geq 98\%$ crystallinity (XRD verified)
- **Passivation** — ASTM A967 nitric acid for stainless steel

Cleanroom Packaging & Sterilisation

Assembly and packaging in ISO 14644 Class 7 cleanrooms. ETO sterilised with biological indicator confirmation. Shelf life validated to 5 years from manufacture date.

Contact

- **Address:** 14 Harley Medical Court, London W1G 9PF, United Kingdom
- **Telephone:** [+91 9319821922](tel:+919319821922)
- **Email:** exports@orthomedimplants.com

This document was generated from <https://www.jksigma.in/manufacturing/>